



## Disinfectant 1 Wipe

### CLASSIFICATION

- ▲ All-in-one bactericide, \*virucide, tuberculocidal, fungicide, deodorizer

### PRODUCT DESCRIPTION

- ▲ *Disinfectant 1 Wipe* is a 1-step Hospital Disinfectant when the use directions for disinfection are followed\*\* (\*\*a pre-clean step required for \*HIV-1, \*HBV, \*HCV, and *C. auris*), and is bactericidal, \*virucidal, fungicidal, and tuberculocidal in the presence of 5% organic soil.

### FEATURES AND BENEFITS

Feature	Benefit
Effective in just 1 minute <sup>‡</sup> against the broad spectrum of organisms listed on the label	1 minute <sup>‡</sup> kill times allow for quick room turnover.
LAS active system	Gentle on surfaces, provides favorable material compatibility.
One-step multi-purpose cleaner & disinfectant	Simplifies daily cleaning and disinfecting processes.
No PPE required	Applicable for EVS room turn over and Nurse's stations.

<sup>‡</sup> Fungi require a 75-second contact time.

### AREAS OF USE

**Disinfectant 1 Wipe is recommended for** use in hospitals, long term care facilities, ambulatory surgical centers, clinics, CAT lab, dental offices, ambulances, emergency rooms, dialysis clinics, home health care, institutional facilities, laboratories, medical offices, operating rooms, patient rooms, public areas, radiology rooms, waiting rooms, day care centers, cafeterias, commercial kitchens, lodging establishments, athletic facilities, exercise facilities, locker rooms, universities, business and office buildings, correctional facilities, laundry rooms, pharmacies, public rest rooms, schools, veterinary clinics, cruise ships, transportation terminals, public facilities, recreational facilities, and retail businesses

#### Use Disinfectant 1 Wipe On:

- ▲ Washable hard, nonporous surfaces of bathroom fixtures, countertops, high touch surfaces, patient monitoring equipment, physical therapy equipment surfaces, tables, bedside tables, computer keyboards, cords, desks, equipment surfaces, faucets, IV stands, light switches, nurse call button, over-bed tables, toys\*\* (\*\*rinse toys with clean, potable water after disinfecting), telephones, exterior toilet bowl surfaces, and tiles, glazed.
- ▲ Non-critical medical device surfaces: wheelchairs, hospital beds, mobile workstations, nonporous mattresses, ultra-sound transducers and equipment, and non-critical equipment surfaces.

- ▲ Other hard, nonporous surfaces made of: acrylic, finished woodwork, Formica®, glazed porcelain, metal, plastic, sealed granite, sealed marble, sealed quartz, sealed stone, stainless steel, vinyl and plastic upholstery. Note: do not use on unpainted wood, natural marble or brass.
- ▲ Use Disinfectant 1 Wipe to clean and disinfect hard nonporous surfaces in health clubs and locker rooms.

#### Disinfectant 1 Wipe:

- ▲ Is formulated for use in daily maintenance programs with a balance of detergents, biocides and malodor counteractants that deliver effective cleaning, disinfecting and malodor control.
- ▲ Will disinfect, clean and deodorize surfaces in rest rooms and toilet areas, sinks and counters, garbage cans and garbage storage areas, and other places where bacterial growth can cause malodors.

Disinfectant 1 Wipe is a 1-step Hospital Disinfectant when the use directions for disinfection are followed\*\*, and is bactericidal, \*virucidal, and tuberculocidal in the presence of 5% organic soil.

\*\* a pre-clean step is required for \*HIV-1, \*HBV, \*HCV, and *C. auris*.

#### **BACTERIA – 1 minute contact time**

*Pseudomonas aeruginosa* ATCC 15442

*Staphylococcus aureus* ATCC 6538

*Bordetella pertussis* ATCC 12743

*Escherichia coli* ATCC 11229

*Escherichia coli* O157:H7 ATCC 35150

*Salmonella enterica subspecies enterica serovar* Typhi ATCC 6539

*Stenotrophomonas maltophilia* ATCC 25596

*Streptococcus pyogenes* ATCC 19615

*Vibrio cholerae* ATCC 14035

#### **ANTIBIOTIC-RESISTANT STRAINS OF BACTERIA – 1 minute contact time**

*Acinetobacter baumannii* (MDR) ATCC BAA-1605

*Escherichia coli* (ESBL) ATCC BAA-196

*Escherichia coli* (NDM-1) CDC 1001728

*Enterococcus faecalis* (VRE) ATCC 51299

*Enterococcus faecium* (MDR) ATCC 51559

*Klebsiella pneumoniae* (KPC) ATCC BAA-1705

*Klebsiella pneumoniae* (MDR) ATCC 51503

*Klebsiella pneumoniae* (NDM-1) ATCC BAA-2146

*Proteus mirabilis* (ESBL) ATCC BAA-856

*Staphylococcus aureus* MRSA ATCC 33592

*Staphylococcus aureus* VISA CDC HIP 5836

*Staphylococcus aureus* CA-MRSA Genotype USA 400 ATCC BAA-1683

#### **VIRUSES\* – 1 minute contact time**

\*Adenovirus Type 5 strain Adenoid 75 ATCC VR-5

\*Cytomegalovirus Strain AD-169 ATCC VR-538

\*Hepatitis B Virus HBV (Duck Hepatitis B as the surrogate)

\*Hepatitis C Virus HCV (Bovine Viral Diarrhea virus as the surrogate) ATCC VR-1422

\*Herpes simplex virus type 1, Strain F(1) ATCC VR-733

\*Herpes simplex virus type 2, Strain G ATCC VR-734

\*Human Coronavirus strain 229E ATCC VR-740

- \*Human Immunodeficiency Virus type 1 \*HIV-1 Strain HTLV-IIIB
- \*Human Rotavirus Strain WA ATCC VR-2018
- \*Human Rhinovirus Type 37 strain 151-1 ATCC VR-1607
- \*Influenza A virus Strain A/Hong Kong/8/68 ATCC VR-544
- \*Influenza B strain B/Lee/40 ATCC VR-1535
- \*Norovirus Feline Calicivirus FCV, strain F-9 as the surrogate ATCC VR-782
- \*Respiratory syncytial virus (RSV) Strain Long ATCC VR-26
- \*Vaccinia virus, Strain WR ATCC VR-119

#### **VETERINARY VIRUSES – 1 minute contact time**

- \*Avian Influenza A strain VNH5N1-PR8/CDC-RG (CDC #2006719965)
- \*Bovine Viral Diarrhea virus (ATCC VR-1422)
- \*Canine Parvovirus Strain Cornell-780916-80 (ATCC VR-2017)
- \*Duck Hepatitis B
- \*Feline Calicivirus strain F-9 (ATCC VR-782)

#### **YEAST – 1 minute contact time**

*Candida albicans* ATCC 10231

#### **FUNGI – 75 second contact time**

*Trichophyton interdigitale* ATCC 9533

*Candida auris* AR-BANK #0381

#### **MYCOBACTERIUM - 45 second contact time @ 22.4°C**

*Mycobacterium bovis* (BCG) (TB) (ATCC 35743)

\*\* a pre-clean step required for \*HIV-1, \*HBV, \*HCV, and *C.auris*.

Disinfectant 1 Spray and Disinfectant 1 Wipe have demonstrated effectiveness against viruses similar to SARS-CoV-2 on hard surfaces. Therefore, these products can be used against SARS-CoV-2 when used in accordance with the directions for use against \*Norovirus and \*Canine Parvovirus at a 1 minute contact time on hard surfaces. Refer to the CDC website at [cdc.gov/coronavirus](https://www.cdc.gov/coronavirus) for additional information.

#### **TO OPEN PACKAGE**

- ▲ Remove cover and discard seal.
- ▲ From center of roll, pull up a wipe corner, twist it into a point, and thread through the hole in cover.
- ▲ Pull through about one inch.
- ▲ Replace the cover.
- ▲ Pull out first wipe and tear off at an angle. (Remaining wipe feed automatically, ready for next use.
- ▲ When not in use, keep cap closed to prevent moisture loss.

#### **DIRECTIONS FOR USE**

- ▲ It is a violation of Federal law to use this product in a manner inconsistent with its labeling.
- ▲ Do not use as a diaper wipe or for personal cleansing.
- ▲ This is not a baby wipe!
- ▲ Do not flush in toilet.
- ▲ For plastic and painted surfaces, spot test on an inconspicuous area before use.
- ▲ Do not reuse or refill empty container.

**DIRECTIONS FOR USE – CONTINUED****TO CLEAN [AND] [&] DEODORIZE HARD, NON-POROUS SURFACES:**

- ▲ Wipe surface until completely wet. Let air dry. Discard wipe in trash.

**TO PRE-CLEAN HARD, NON-POROUS INSTRUMENTS PRIOR TO TERMINAL STERILIZATION/HIGH-LEVEL DISINFECTION:**

- ▲ Instruments must be thoroughly pre-cleaned to remove excess organic debris, rinsed, and rough dried. Clean and rinse lumens of hollow instruments before wiping with this product. Wipe all surfaces of instruments with this product until thoroughly wet. Allow to remain wet for listed contact time. Wipe with a clean, damp cloth or paper towel and allow to air dry.

**TO CLEAN AND DISINFECT HARD, NON-POROUS SURFACES IN A VETERINARY APPLICATION:**

- ▲ For disinfection of pre-cleaned animal quarters and kennels. Remove all animals and feed from premises, animal transportation vehicles, crates, etc. Remove all litter, droppings, and manure from walls, floors, and surfaces of facilities occupied or traversed by animals. Empty all feeding and watering equipment. Pre-clean all surfaces with soap or detergent, and rinse with water. Wipe surfaces and let stand for required contact times. Ventilate buildings and other closed spaces. Do not house animals or employ equipment until treated surfaces have been thoroughly rinsed with water and allowed to dry. Thoroughly scrub all treating, feeding, and watering appliances with soap or detergent, and rinse with potable water before re-use.

**TO CLEAN AND DISINFECT HARD, NON-POROUS SURFACES:**

- ▲ Wipe surface until completely wet. Allow the treated surface to remain wet for stated contact time. Let air dry. For visibly soiled areas, a pre-cleaning is required.

**KILLS \*HIV-1, \*HBV, AND \*HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS**

- ▲ This product kills \*HIV-1, \*HBV, and \*HCV on pre-cleaned environmental surfaces/objects previously soiled with blood/body fluids in health care settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of \*HIV-1 \*Hepatitis B virus or \*Hepatitis C virus.

**SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST \*HIV-1, \*HBV, AND \*HCV OF SURFACES/ OBJECTS SOILED WITH BLOOD/BODY FLUIDS**

- ▲ PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable impervious gloves, gowns, masks, and eye coverings.
- ▲ CLEANING PROCEDURE: Blood and other body fluids must be thoroughly cleaned from surfaces and other objects before application of this product. This cleaning may be accomplished with any cleaning solution, including this product.
- ▲ CONTACT TIME: Thoroughly wet the surface with Disinfectant 1 Wipe. Allow surface to remain wet for 1 minute. The efficacy of a 1 minute contact time has been shown adequate against \*HIV-1 \*HBV, and \*HCV.
- ▲ DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids must be autoclaved and disposed of in accordance with federal, state, or local regulations.

**SPECIAL LABEL INSTRUCTIONS FOR CLEANING PRIOR TO DISINFECTION AGAINST CANDIDA AURIS:**

- ▲ **PERSONAL PROTECTION:** Wear appropriate barrier protection such as gloves, gowns, masks, or eye covering.
- ▲ **CLEANING PROCEDURE:** Fecal matter/waste must be thoroughly cleaned from surfaces/objects before disinfection by application with a clean cloth, mop, and/or sponge saturated with the product. Pre-cleaning is to include vigorous wiping and/or scrubbing and all visible soil is removed. Surfaces in patient rooms are to be cleaned in an appropriate manner, such as from right to left or left to right, on horizontal surfaces, and top to bottom, on vertical surfaces, to minimize spreading the organism. Restrooms are to be cleaned last. Do not reuse soiled cloths.
- ▲ **INFECTIOUS WASTE DISPOSAL:** Materials used in the cleaning process that may contain feces/wastes are to be disposed of immediately in accordance with local regulations for infectious materials disposal.
- ▲ This product is not to be used as a terminal sterilant/high-level disinfectant on any surface or instrument that is (1) introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean critical or semi-critical medical devices prior to sterilization or high-level disinfection.

## SPECIFICATIONS

<b>Ingredients</b>	Dodecylbenzene sulfonic acid	0.29%
	Other Ingredients:	99.71%
	Total:	100.00%
<b>Appearance</b>	Absorbed liquid (wipes)	
<b>Color</b>	White (liquid: clear, colorless)	
<b>Odor</b>	Sweet	
<b>pH</b>	2.1 – 2.4	
<b>Flammability</b>	Not flammable, does not sustain combustion	
<b>Relative density</b>	1.007 – 1.009	
<b>Product Expiration</b>	365 days	

## ENVIRONMENT

- ▲ Hazard Rating: Health – 0; Flammability – 0; Physical Hazard – 0

## REGULATORY REFERENCES

- ▲ US EPA Reg. No. 1677-263

**STORAGE AND DISPOSAL**

- ▲ Store in original container in areas inaccessible to children and persons unfamiliar with its proper use.
- ▲ **PRODUCT STORAGE:** Product must be stored in a dry area, away from direct sunlight and heat. Do not store in direct sunlight.
- ▲ **TOWELETTE DISPOSAL:** Do not reuse towelette. Dispose of used towelette in trash. Do not flush in toilet.
- ▲ **CANISTER:** Non-refillable container. Do not reuse or refill empty container. Rinse and recycle empty container or discard in trash.
- ▲ **KEEP FROM FREEZING**

**PRECAUTIONARY STATEMENTS****HAZARDS TO HUMANS AND DOMESTIC ANIMALS**

**CAUTION:** Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet.

**FIRST AID**

**IF IN EYES:** Hold eye open and rinse slowly and gently with water for 15-20 minutes. Remove contact lenses, if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

**FOR EMERGENCY MEDICAL INFORMATION CALL TOLL FREE: 1-800-328-0026  
OUTSIDE NORTH AMERICA, CALL 1-651-222-5352**

Have product container or label with you when calling a poison control center, doctor, or going for treatment.

**PACKAGING & PRODUCT NUMBERS**

EPA Reg. No. 1677-263

Size Available	Product Number
220 ct x 12 cans (wipe size: 6" x 6.75")	6000261
85 ct x 12 cans (wipe size: 9" x 12")	6000262

<b>Worldwide Headquarters</b> Ecolab   1 Ecolab Place   St. Paul, MN 55102 USA <a href="http://www.ecolab.com">www.ecolab.com</a> 1.866.781.8787	
---	--