
OPTIPRO™ GEL INSTRUMENT PRE-CLEANER

INSTRUMENT REPROCESSING



A ready-to-use surgical pre-cleaner that helps remove heavy soils from surgical instruments, maintains moisture during transport, and protects valuable instruments.

Strong Performance

- Advanced gel technology clings to heavy soils and keeps them from drying on instruments
- Dissolves blood and fat, making removal from instruments fast and effective
- ▲ Neutral pH allows for safe use with most instrumentation material

Superior Technology

- ✓ Non-aerosol pressurized Bag-on-Valve container uses inert compressed air to dispense product*
- Blue color makes it easy to ensure uniform coverage of instruments
- Gel form avoids splashing of liquid during transport
- ▲ Ready-to-use, convenient and portable for OR application

Sustainable Advantage

▲ Phosphate-free formula to help reduce environmental impact

OptiPro™ Gel Instrument Pre-Cleaner



SKUs

Size	Description	Part #
12 - 32 oz.	Trigger Spray	6075552
6 - 13.7 oz.	Bag-on-Valve	6000252



^{*}No propellant mixes with product