



THE QUICK AND EASY ALTERNATIVE TO MOP AND BUCKET SYSTEMS



Klerwipe™ Neutral Detergent Mop Wipes

Non-woven, low particulate wipes impregnated with a neutral, non-ionic detergent

Ecolab's Klerwipe Neutral Detergent mop wipes are a quick and convenient method for cleaning large surfaces. They are designed for use with the Klercide autoclavable mop wipe frames for optimum surface coverage and remove the need to use space-consuming buckets or trolley systems.

Klerwipe Neutral Detergent Mop Wipes are:

- ▲ Hydroentangled polyester/cellulose
- ▲ Impregnated with Klercide Neutral Detergent
- ▲ Non-woven and binder free
- ▲ Basis weight: 58-78gsm

PRODUCT QUALITY ASSURANCE

Manufactured in our state-of-the-art facilities to the principles of cGMP, Klerwipe Neutral Detergent pouch wipes are safe for use in core manufacturing areas.



Klerwipe™ Neutral Detergent Mop Wipes



Non-woven, low particulate wipes impregnated with a neutral, non-ionic detergent

Why use Klerwipe Neutral Detergent Mop Wipes?

EASE OF USE

Ready-to-use, quick and easy to apply, Klerwipe Neutral Detergent Mop Wipes are impregnated with Klercide Neutral Detergent, a low-foaming, non-ionic solution with neutral pH. These non-corrosive wipes are powerful and effective in a wide range of cleanroom applications.

Klerwipe Neutral Detergent Mop Wipes are designed to fit easily onto the Klercide Mop Wipe Frames. Compared to bucket systems, they are:

- ▲ Faster to use
- ▲ Easier to manage
- ▲ Able to eliminate the potential risk of a dilution error when mixing a concentrate

CONSISTENCY

Ecolab Klerwipes are highly absorbent and remain wet throughout use, ensuring a consistent application of detergent to the entire wiping surface area. The use of wiping can increase cleaning consistency and the quality of the final product.

PRODUCT QUALITY

Klerwipes are manufactured in Ecolab's ISO 9001 accredited facilities to the principles of cGMP. These manufacturing conditions ensure Klerwipe Neutral Detergent wipes are manufactured to the highest quality and are compliant for use in your cleanroom areas.

The Klerwipes are:

- ▲ Impregnated with Klercide Neutral Detergent
- ▲ Made with Water for Injection (WFI)
- ▲ Packed and bagged in an ISO Class 7 cleanroom
- ▲ Supported with Certificates of Analysis and Certificates of Sterility with every batch
- ▲ Supported with full documentation



Klerwipe Neutral Detergent
420 x 250mm

PHARMACEUTICAL

Safety. Compliance. Efficiency.

APPLICATION AREA

For use on floors, walls and other large, hard surfaces.

INSTRUCTIONS FOR USE

For use with the Ecolab Klercide mop wipe frames. Pull back the upper layer of the peel-pouch pack and place wipe frame head on the mop wipes, wrap the wipe around the head and secure either by pushing firmly into plugs, or with stainless steel clips.

INGREDIENTS

Fatty alcohol ethoxylate, tetrahydropropyl ethylenediamine and Water for Injection.

PACKAGING INFORMATION

Neutral Detergent Mop Wipes

SUBSTRATE	UNIT SIZE	UNIT OF SALE
45% Polyester/ 55% Cellulose	420 x 250mm mop wipe	10 x 20 wipes

FOR PRODUCT ORDER CODES AND ADDITIONAL INFORMATION, PLEASE CONTACT YOUR ACCOUNT MANAGER OR VISIT OUR WEBSITE

SHELF LIFE

Unopened, this product has a one year shelf life.

ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences
www.ecolab.com/pharmaceutical