

P3-cosa FLUX

Cleaning of membrane systems
in the Pharmaceutical, Cosmetic and
Biotechnological Industry



- Adapted for membrane cleaning validation
- All products are both liquid and water-soluble
- Produced under GMP

P3-cosa FLUX range

Product	Short Description	
P3-cosa FLUX 22	neutral	surfactant-based cleaner, to be used alone or in combination with neutral, alkaline or acidic products from the P3-cosa FLUX range
P3-cosa FLUX 33	alkaline	strong alkaline, surfactant-free cleaner, to be used alone or in combination with P3-cosa FLUX 22 or P3-cosa FLUX 48
P3-cosa FLUX 44	alkaline	mild alkaline, surfactant-free cleaner, to be used alone or in combination with P3-cosa FLUX 22 and/or P3-cosa FLUX 48
P3-cosa FLUX 48	neutral	enzymatic, surfactant-based cleaner, to be used alone or in combination with neutral or alkaline products from the P3-cosa FLUX range
P3-cosa FLUX 55	acidic	strong acidic, surfactant-free cleaner for the removal of any mineral scaling of water hardness and phosphates, to be used alone or in combination with P3-cosa FLUX 22

Detailed product recommendations are available via our CD-ROM "P3-Membrane Care" or can be given on special request.

In order to help us better understand your requirements, please consider the following when requesting additional information:

Membrane type: spirral, tube, flat sheet ...
 Material: polysulphone, ceramics ...
 Limits: pH, chlorine, temperature ...
 Product: fermentation broth, antibiotics, blood ...
 Problem: capacity, retention ...

