


MULTI-USE SERIES 3003 CIP UNIT

Product Description

The Ecolab Model 3003 CIP unit is skid mounted, pre-piped, pre-wired and ready for owners hook-up of steam, water, air, electricity and piping.

Benefits

- Product Quality
- Conserves Water
- Program Flexibility
- Personnel Safety
- Multi-Functional
- Cost Efficient
- Chemical Control



Standard Features

- ISO 9002 Manufacturing
- Meets 3-A Std. #605-04
- Fristam Supply Pump
- Waukesha Valves
- ASCO Utility Valves
- Supply Line Strainer
- Supply Line Check Valve
- Return Check Probe
- Schedule 10, 304 SS Frame
- 120 Gallon, 12 Gauge Water Tank
- Foxboro Flow Meter
- 240 Gallon, 12 Gauge Recovery Tank

- 240 Gallon Caustic Tank
- Enerquip SS U-Tube Heat Exchanger
- Foxboro Conductivity Transmitter
- Foxboro Pressure Transmitter
- Manifold Recirculation Valve
- 3-Level Probes
- 3-Point Chemical Injection
- Sample Valve

Applications

The 3003 CIP System is specified for dairy plants where water conservation and operational costs are of primary importance, and program flexibility is required for multiple cleaning applications.

