


PocketSwab[®] Plus

Shelf Stable



The PocketSwab[®] Plus is a self-contained single-service, ATP rapid hygiene test. Food processors recognize the importance of cleaning and sanitizing as part of their Quality, SSOP and HACCP programs. PocketSwab Plus enhances food safety programs as part of a total brand protection strategy, with immediate feedback to monitor and guide the sanitation process.

Improved Shelf Stability

- No refrigeration required.
- Room temperature stable for 12 months.

Superior Sensitivity

- Up to 40 times more sensitive than other ATP-based rapid tests.
- * Manufactured under strict clean room conditions.

Consistent, Reliable Results

- Proven repeatability — laboratory testing in swab variability under standard, controlled conditions.
- ATP positive control available to check system calibration and swab integrity.

Quick and Easy to Use

- Real-time results on the plant floor.
- Metal detectable, write-on labels.

Unique, Patented Design

- Thread design assures consistent results from user to user.
- Unique luciferin / luciferase tablet dissolves on contact for improved accuracy and stability.
- Clear microtube allows operator to observe reaction.
- PocketSwab Plus contains no glass or hazardous substances — may be discarded in ordinary waste containers.

PocketSwab Plus

Shelf Stable

Monitoring of Sanitation Effectiveness
has never been Easier or Faster