

# NOVA LUM™



*All novaLUM tests are conveniently stored, tracked and trended by the new novaLINK™ software.*

The novaLUM™ is a revolutionary, lightweight, ergonomically designed luminometer to aid in ATP-based hygiene monitoring programs in the dairy, food and beverage processing industries. A wireless option is available to transmit data in real time.

## **Sensitive**

- Flexible for low, medium and high risk processing environments
- Only system with ATP indicator for allergen cross contamination presence

## **Sleek**

- Ergonomic design with direct swab entry
- One hand portability
- Smallest, lightest luminometer with photomultiplier tube

## **Smart**

- Intuitive yet powerful novaLINK™ software
- Flexible reporting function allows customized reporting

## **Streamlined**

- Fastest pre-op hygiene test available — 5 second read time
- Full numeric keypad
- Rocker scroll keys for quick navigation

## **Secure**

- Track and trace up to 6000 test results
- Manage multiple sampling plans — 400 test sites per plan

**ECOLAB®**