



MODEL 4412 DROP STATION LOW PRESSURE - CENTRAL FOAM

Product Description

The Ecolab Model 4412 Drop Station is a wall mounted low-pressure, central foam station, ready for owner's hook-up of central detergent supply and air. It generates foam by mixing compressed air into a detergent stream that is supplied from a central detergent dispensing system.

Benefits

- Compact
- Durable
- Easy to Use
- Cost Efficient
- Controlled Detergent Concentrations

Standard Features

- Foam injector
- SS construction
- Mounting brackets
- Shut-off ball valves
- Air Regulator
- Quick connect hose fitting
- SS mixing chamber
- SS heavy duty hose hanger

Applications

The 4412 low pressure central foam drop station is for generation of detergent cleaning foam using detergent supplied from a central dispensing system, and mixing it with compressed air injection. This provides convenient manual foam application for facility and product contact surfaces in food or beverage processing facilities.

