



# MODEL 4414 DROP STATION LOW PRESSURE - CENTRAL SANITIZE

## Product Description

The Ecolab Model 4414 Drop Station is a wall mounted low-pressure, central sanitize station, ready for owner's hook-up of central sanitizer supply. It delivers sanitizer that is supplied from a central sanitizer dispensing system.

## Benefits

- Compact
- Durable
- Easy to Use
- Cost Efficient
- Controlled Sanitizer Concentrations

## Standard Features

- SS construction
- Mounting brackets
- Shut-off ball valve
- Quick connect hose fitting
- SS heavy duty hose hanger

## Applications

The 4414 low pressure central sanitizer drop station is for dispensing of sanitizer at the proper concentration, supplied from a central sanitizer dispensing system. It will provide total coverage of equipment surfaces without atomizing the sanitizer into the air for convenient manual sanitizer application on facility and product contact surfaces in food or beverage processing facilities.

