

SPRAY-IT PORTABLE SPRAY UNITS 5 AND 20 GALLON



5 GALLON SPRAY UNIT



20 GALLON SPRAY UNIT

Primary Applications:

Portable units used for sanitizing equipment, walls, floors, and other environmental surfaces

Features:

- Draws from pre-mixed solution
- 30 ft. (9 m) of discharge hose, with 3/8 in. (9.5 mm) inside diameter
- 10-12 ft. (3.0-3.7 m) output distance
- 2 gal/min (7.6 l/min) flow rate
- Polyethylene tank with hinged back plate
- Natural color tank for visibility of product levels
- Santoprene™ diaphragm pump

Utilities Required:

- Compressed Air: 40-80 psi (3-5 bar) with 2 cfm (28.3 l/min)
- Liquid Temperature Range: 40-100°F (4.4-37°C)

Limitations:

- Limited by the length of the air and chemical hoses

Dimensions/Weight:

5 Gallon: 24" x 15" x 21" and 28 lbs (empty)

20 Gallon: 22" x 17" x 47" and 42 lbs (empty)

What's in the box?

- Complete foam unit with discharge hose and wand
- I/O manual

Replacement Parts:

For part numbers, see Repair Parts Manual, 'SANITIZING EQUIPMENT' section.

SPRAY-IT is a brand of Innovative Cleaning Equipment Inc.
Santoprene™ is a registered trademark of Exxon Mobil

Application Equipment → Sanitizing Equipment → Sprayers

Worldwide Headquarters
370 Wabasha Street N St. Paul, MN 55102
www.ecolab.com 1 800 392 3392

© 2015 Ecolab USA Inc. All rights reserved. 43485/0300/0112

