

Helping You Improve Restroom Cleanliness.

Ecolab's time-saving solutions make it easy to offer disinfected restrooms all day, every day.

For every restroom cleaning challenge, Ecolab has a solution.

CHALLENGE	SOLUTION
Time required to keep restrooms clean	The Ecolab Cleaning Caddy provides touch-free cleaning in half the time.
Training staff on proper restroom cleaning procedures	Ecolab provides effective, easy-to-use solutions and training tools to help get the job done quickly and easily.
Managing product inventory and reducing overall cost	Ecolab's concentrated cleaners and innovative dispensing systems optimize product use and minimize waste.



The Cleaning Caddy and QC™ products work together to deliver daily disinfecting with superior ease-of-use.

- ▲ Closed-loop dispensing to avoid free pouring
- ▲ Concentrated with less packaging waste to reduce environmental impact

PRODUCT NAME	DESCRIPTION	PACK SIZE	PICK CODE
MULTI-PURPOSE DISINFECTANTS			
 QC 57 PEROXIDE MULTI-SURFACE CLEANER & DISINFECTANT	A 2 in 1 product that is a hard surface cleaner and glass cleaner. Peroxide based, disinfectant with fast contact times (U.S. between 5 and 3 minutes, Canada 5 minutes). Effective against Norovirus. U.S. EPA Reg. No. 1677-238. Canada DIN 02412314.	2x1.3L	6100792 (US Only) 6100896 (Canada Only)
 NEUTRAL DISINFECTANT CLEANER	Multi-purpose, germicidal detergent that disinfects, cleans and deodorizes in one labor-saving step. Effective on hard surfaces and can also be used on non-porous athletic mats, exercise and personal safety equipment. U.S. EPA Reg. No. 47374-129-1677. Canada DIN 02342820.	2x1.3L	14541 (US Only) 6100389 (Canada Only)
 A-456® II DISINFECTANT CLEANER	Broad spectrum, one-step, hard surface, disinfectant, cleaner and sanitizer. Delivers effective cleaning and biocidal action against bacteria, fungi and viruses. Effective against Norovirus. EPA Reg. No. 6836-78-1677.	2x1.3L	61166931 (US Only)
MULTI-PURPOSE CLEANER			
 QC 54 PEROXIDE GLASS & SURFACE CLEANER	Daily multi-purpose cleaner gives a streak-free shine through the use of peroxide. Formulated to clean windows, countertops, glazed porcelain, metal fixtures, walls, impervious floors and other hard surfaces.	2x1.3L	6100610 
GLASS CLEANER			
 QC 50 ULTRA CONCENTRATED GLASS CLEANER	Streak-free, fast-drying performance on windows, mirrors, glass, plastic, acrylics and painted surfaces. Helps remove tough soils including smoke and grease films. Low sudsing.	2x1.3L	6100291 
BATHROOM CLEANERS			
 73 DISINFECTING ACID BATHROOM CLEANER	Acid disinfecting cleaner helps remove hard water build-up and soap scum from tubs, showers, countertops, toilets and sinks. Disinfects with reduced acid flashback. EPA Reg. No. 1677-204.	2x1.3L	6101039 (US Only)
 QC 92 ULTRA CONCENTRATED NEUTRAL BATHROOM CLEANER	Neutral cleaner formulated to help remove soils from toilet bowls, urinals, and shower walls and floors.	2x1.3L	6100280 
 QC 93 ULTRA CONCENTRATED ACID BATHROOM CLEANER	Acid bathroom cleaner helps effectively remove dirt, soap scum and hard water deposits. Safe for use on bathroom chrome, porcelain, vinyl, tile and fiberglass.	2x1.3L	6100295 
ODOR CONTROL			
 BIO-ENZYMATIC ODOR ELIMINATOR (QC 77)	An all-in-one product that instantly knocks down odor upon application and has continuous action to eliminate odor at the source. It is specifically designed to counteract urine and body odor, effective on hard and soft surfaces with a pleasant fragrance.	2x1.3L	6101023
FLOOR CLEANERS			
 QC 34 HIGH PERFORMANCE ULTRA CONCENTRATED NEUTRAL FLOOR CLEANER	Neutral floor cleaner designed for daily use in automatic scrubbers and mop buckets. Powerful blend of surfactants help remove soils without leaving a film or haze. Cleans without removing finish or dulling gloss and does not require rinsing. Safe on most floor types. Very low sudsing.	2x1.3L	6100082 
 QC 35 MAXX DUAL ACTION FLOOR CLEANER (BRIGHT SPEED)	Dual-purpose alkaline floor cleaner for daily maintenance and scrub and recoat use. Clean, fresh scent is less disruptive for daily cleaning while ultra-low foaming technology increases auto scrubber productivity. Very low sudsing.	2x1.3L	6100433



ECOLAB CLEANING CADDY (QC™)

Rethink Restroom Cleaning

Daily Disinfecting • Better Results • Lower Total Cost

ECOLAB®
Everywhere It Matters.™

370 Wabasha Street N St. Paul, MN 55102
www.ecolab.com 1 800 35 CLEAN

©2016 Ecolab USA Inc. All rights reserved. S45508/0402/0616 9265-2509

ECOLAB®
Everywhere It Matters.™

Better Results, Lower Costs.

Improve the cleaning process using a better cleaning system and enjoy the results.

Cleaner Restrooms

- Clean and odor-free restrooms create a great first impression
- Touch-free system simplifies the cleaning process

ECOLAB CLEANING CADDY GETS RESTROOMS
54% CLEANER
THAN TRADITIONAL CLEANING METHODS*

Everyday Disinfecting

- Effective disinfectant solutions help prevent the spread of germs
- Deep clean and disinfect every day of the week vs deep cleaning weekly/monthly

Operational Efficiency

- Increase employee productivity by cleaning in half the time**
- Consistent cleaning process for better results



Ecolab Cleaning Caddy

Touch-Free Restroom Cleaning

Smart Design

- Closed-loop product dispensing
- Spray wand (90 PSI spray pressure) evenly distributes product
- Rechargeable battery lasts 2.5 hours on a single charge
- On-board tool storage
- Wet vacuum removes dirty solution for a complete clean**

Easy to Use

- Easy-access tank holds up to 13.5 gallons of water
- Compact and easy to maneuver, similar footprint to mop and bucket
- Easy-to-use control panel simplifies operation



Cleaning Caddy Control Panel

Caddy + Wet Vacuum

9265-1053

Includes:

- On-Board Wet Vacuum
- Glass Squeegee
- Hand Caddy
- Doorstops (2)
- Dual Surface Deck Brush
- Charger
- Welcome Packet & Training Material



Caddy + Shelf

9265-1039

Includes:

- Storage Shelf
- Glass Squeegee
- Hand Caddy
- Doorstops (2)
- Dual Surface Deck Brush
- Charger
- Welcome Packet & Training Material



Cleaning Caddy Accessories

ITEM NAME	DESCRIPTION	PICK CODE
Wet Vacuum		
WET VACUUM	10-gallon wet vacuum kit	9265-2306
WET VACUUM HEAD	Wet vacuum replacement head	9265-2997
WET VACUUM SQUEEGEE BLADE REPLACEMENT KIT	Wet vacuum squeegee blade replacement kit	9265-2322
WET VACUUM FILTER REPLACEMENT	Wet vacuum filter replacements (3-pack)	9265-2317
Cleaning Tools		
BROOM HANDLE-BLUE	Combo scrub brush & squeegee handle	60302-01-00*
COMBINATION SCRUB BRUSH/SQUEEGEE HEAD	Combo tool scrub brush & squeegee head	60334-01-00*
DUAL-SURFACE DECK BRUSH-BLUE	Replacement deck brush head	8999-0051
FIBERGLASS INTERLOCKING HANDLE-BLUE	Dual-surface deck brush handle	8999-0045
TOILET BRUSH WITH HOLDER	Toilet brush tool and holder (12-pack) Replacement Mops (12-pack)	9263-1324 9201-1186
GLASS SQUEEGEE	12" glass squeegee (3-pack)	9621-1112
GLASS SQUEEGEE REPLACEMENT BLADE	Glass squeegee replacement blade	60332-10-00*
WORK SURFACE SQUEEGEE-BLUE	6.5" handle, 9" blade	30413-01-00*
Miscellaneous		
"DISINFECTED DAILY" STICKER-WHITE	White "Restroom Disinfected Daily" sticker (1 sheet of 6)	9265-2329
"DISINFECTED DAILY" STICKER-CLEAR	Clear "Restroom Disinfected Daily" sticker (1 sheet of 6)	9265-2328
DOORSTOP	Doorstop	9265-2327
HAND CADDY	Hand caddy holds tools and up to 3 spray bottles	8999-2948
MOP SINK HOSE	6' hose	9480-2188

* Results based on over 1,000 swab samples collected and analyzed by Ecolab.

** Based on ISSA cleaning times for restroom floor cleaning.

** Power cord is required to operate wet vacuum.