


# Quality Water MAKES A DIFFERENCE



## INTRODUCING THE KAY® QSR MANAGED WATER FILTRATION PROGRAM

Grounded in cutting-edge technology, backed by Kay's trusted expertise  
and delivering the best beverage experience for your guests

KAY® PROGRAM

**ECOLAB®**

# Quality Water EXPERTLY DELIVERED

---

A comprehensive, managed water filtration program that puts the focus on beverage quality

---

## INDUSTRY LEADING TECHNOLOGY

- Patented design for high flow, high capacity applications
- Patented hollow carbon media
- Chlorine and chloramines removal
- Cyst reduction

## PROVEN SERVICE AND EXPERTISE

- Customizable platform based on beverage applications and local water conditions
- Automatic annual filter change outs completed by Kay Territory Representative
- Unit level water testing and reporting to ensure consistency and quality
- Filtration system inspections during Value Visits





# Quality Water MATTERS

**71%** OF CONSUMERS NOW PURCHASE  
A BEVERAGE AWAY FROM HOME TWICE A WEEK OR MORE\*

The KAY® QSR Managed Water Filtration Program aligns **cutting edge technology** with industry leading **service and expertise** – a complete solution that is critical to **protecting the quality** of the beverages that you serve.

## ENHANCE GUEST SATISFACTION

- Deliver consistent, high quality beverages
- Reduce poor tastes and unpleasant odors
- Improve clarity in teas and other beverages

## IMPROVE OPERATIONAL EFFICIENCY

- Reduce energy consumption
- Protect against expensive service calls and equipment downtime
- Optimize equipment life and performance

## SUPPORT FOOD SAFETY

- Sanitary, quick-change technology
- Protection against parasitic cysts



*The KAY QSR Managed Water Filtration Program includes everything from installation to filter change outs – allowing you to focus on what matters most—your guests.*

\*Technomic's "Beverage Consumer Trend Report," 2012



# Customizable technology to meet your needs

- Automatic annual change outs
- High flow, high capacity
- Standard 0.5 micron particulate reduction
- Sanitary, quick change technology
- Removes both chlorine and chloramines



## FILTER CARTRIDGES

SIZE		10"	14"	20"
SPECIFICATIONS - CHLORINE	FLOW (GPM)	5.0	7.0	10.0
	CAPACITY (GALLONS)	50,000	75,000	100,000
SPECIFICATIONS - CHLORAMINES	FLOW (GPM)	2.0	3.7	3.0
	CAPACITY (GALLONS)	8,000	15,000	22,500
PART NUMBER/MODEL		92213684 KAY-T010	92213686 KAY-T014	92213687 KAY-T020
SCALE CARTRIDGE PART NUMBER/MODEL		92213688 KAY-T010S	92213689 KAY-T014S	92213690 KAY-T020S

## EQUIPMENT SYSTEMS



**SINGLE**  
Part #: 92213705



**TWIN**  
Part #: 92213707



**SPLIT-TWIN**  
Part #: 92213708



**TRIPLE**  
Part #: 92213710

For more information or technical assistance:

Call 800-529-5458 in US or Canada

or +1 336-931-2000 worldwide