

A Hard-Working Ergonomic Laundry Program



Ecolab Solid Laundry Program

- ▲ Enhance guest and residence safety with skin friendly pH and sanitization
- ▲ Advance simplified procedures and training
- ▲ Help promote safety with sanitization
- ▲ Enable consistent results and operational effectiveness with out-of-product alarms

*Ecolab Laundry Program, along with Ecolab laundry experts and a strong focus on stain removal and rewash reduction, can potentially **help reduce your rewash rate to 3-5%***

¹ According to Rewash: Money Down the Drain published in 2017 by American Laundry News, the national average reject rate is 3-5%.



Ecolab Solid Laundry Program



Tackles tough stains, which is the biggest challenge of laundry professionals¹ and the #1 reason for discarded linens²

- High-performing stain fighting removes stains and restores whiteness and brightness across various laundry conditions

Promotes Safe, Simple Operation for Staff

- Requires less physical strain compared to traditional liquid laundry programs
- Closed-loop system reduces exposure to chemicals
- Easy-to-use and simplifies training
- Help reduce rewash and enable consistent results with out-of-product alarm

22% of hotels reported injuries involve improper lifting and carrying heavy objects³

¹ 77% of laundry professionals in the lodging industry said stain removal presents their biggest challenge. 2012 Ecolab Laundry Habits and Practices Study. Conducted by Satisfaction Management Systems. Reported in February 2013.
² Ecolab in Market Linen-Life Research, June 2014.
³ Loss Prevention COMPInsiders, Top Injuries in Hotels and Motels, Summit Consulting LLC, 2020.

SIMPLE 3-PRODUCT PROGRAM

Solid Laundry Detergent

Enables consistent results in various laundry conditions



- Now with improved detergency
- Phosphate and NPE free - helps protect waterways and the environment

Solid Laundry Sour Soft

Helps prevent linen damage and improves user experience



- Now with new fresh fragrance, iron control, and antichlor
- No need for add-on product

Low Temp Laundry Solid Chlorine Sanitizer



- Solid Sanitizer provides excellent bleaching and helps preserve linen integrity and strength
- Sanitizes linen at the right dosage

PRODUCT	SKU	PACK SIZE
Ecolab Solid Laundry Detergent	6102325	4 - 9 lb
Ecolab Solid Laundry Sour Soft	6102327	2 - 6 lb
Ecolab Low Temp Laundry Solid Chlorine Sanitizer	6100975	2 - 4 lb

For more information, **contact your Ecolab Representative, call 1 800 35 CLEAN, or visit Ecolab.com.**