

ECO-FLEX™

THE LATEST INNOVATION IN UDDER CARE

Fast, effective, and affordable formula
designed for teat disinfection before
and after milking

Complete 3 in 1 TECHNOLOGY

- ✓ **KILLS**
- ✓ **PROTECTS**
- ✓ **CONDITIONS**



POWERFUL DISINFECTION

- Kills **>99.999%** of mastitis-causing pathogens, including *S. aureus*, *E. coli*, *Pseudomonas*, *S. dysgalactiae*, *S. agalactiae*, *S. uberis*, *K. pneumoniae*, and *E. aerogenes*



EFFICIENT PROTECTION

- Delivers **>6 log reduction in 15 seconds** including wild strains of *S. aureus*
- Kills *S. aureus* up to **4X faster** than a 1% iodine treatment



HELPS PROMOTE TEAT SKIN HEALTH

- **Pre-milking emollient** delivers **4% conditioning**
- **Post-milking emollient** delivers **8% conditioning** system, a unique blend of glycerin, sorbitol, lanolin, and propylene glycol

AN ALL-NEW, TRIPLE-ACTIVE ANTIMICROBIAL SYSTEM

Pelargonic Acid

- Organic acid, naturally found in plants/animals; commonly used for fruit & vegetable treatment

Lactic Acid

- Organic acid, naturally found in milk and skin; used as skin moisturizer, exfoliant & barrier aid


Dodecylbenzene Sulfonic Acid

- Effective "clean & kill" treatment

ECO-FLEX™

THE LATEST INNOVATION IN UDDER CARE


AOAC Test - *S. aureus* 10% Milk Complete Kill



POWERFUL DISINFECTION

Innovative formulation enables **fast, >6 log reduction kill** with wild strains of *Staphylococcus aureus* under challenged conditions


AOAC Test - *S. aureus* 10% Milk



EFFICIENT PROTECTION

Triple active ingredient system works fast to kill pathogens **4x faster** than iodine treatment

Teat Health: Post-Milking Trials



HELPS PROMOTE TEAT SKIN HEALTH

Unique blend of compounds delivers post-milking conditioning, and may contribute to the **improvement of teat skin health**